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Digitalwaage

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Operating Instructions Electronic Hanging Balance

E

KERN HCB	20K50	50K100	100K200	200K500
Read-out	50 g	100 g	200 g	500 g
Weighing Range	20 kg	50 kg	100 kg	200 kg

To fit or exchange the battery

Push off the battery cover on the back of the balance in the direction of the arrow. Insert **3 AAA** batteries. Replace the battery cover. When the battery power is used the display will show "LO". Press the **OFF** key and change the batteries at once. To preserve battery life the balance switches off automatically after **4** minutes if no weighing has taken place.

Weighing

- The balance is switched on by pressing the "ON/OFF" key.
- The balance display shows "0.00".
- Hook the object being weighed.
- The steady reading is the weight of the object.
- If the weighing object is heavier than the weighing range the display indicates "E" (=Overload).
- The balance is switched off by pressing the "ON/OFF" key once.

Taring


- Switch the balance on with "ON/OFF" key and wait for "0" to appear in the display. Hang the tare load and press the "TARE" key again.
- The balance display becomes "0.00". The weight of the tare load has however been memorised internally. The item to be weighed can now be placed in the tare container and the weight can be read off..

Data HOLD

- When the **HOLD** key is pressed, the display value (even when the weight values vary) shown on the display will be held 5 seconds.

Important using informations

NO


- Don't use with pressure

NO

- Don't twist

O.K.

- Just using with tension

Important Notices

These electronic scale is a precision instrument. Electromagnetic fields can cause major display discrepancies. The scales must then be repositioned away from electromagnetic fields. All sources of environmental interference, such as drafts and vibrations, should be avoided. Sudden changes of temperature should be avoided. The scales must be reset to match changes in temperature.

The scale is not hermetically sealed, therefore avoid high humidity, steam and dust. Do not bring liquids into direct contact with the scales, as these can penetrate into the measuring mechanism. Cleaning material should only be dry or barely damp. Do not use solvents as these can damage paintwork or other plastic parts. Remove damaged items immediately from the scale.

The measuring mechanism will be stabilised by allowing the scale to warm up for a few minutes after switching them on. Place items to be weighed carefully on the scales.

Do not hang objects under the scale for any period of time, apart from normal use. Sudden shocks or overloading the scales beyond the maximum permitted weight should absolutely be avoided, scale could be damaged.

In case of problems operating the weighing program, switch the scales briefly off and on. The weighing set-up must then be restarted right from the beginning.

Never operate the scales in areas where there is a danger of explosion, the models in this series are not protected against explosion.

Check the scales regularly against, known, external test weights.

Opening the scales or failing to use them in accordance with the written instructions will invalidate the warranty.

Please keep all packaging material for possible return of the scales. Scales must only be returned in their original packaging.

Safety warning

The hanging electronic scales from the KERN HCB series are only suitable for hanging from a mechanical hook, e.g. a crane hook. Never use the balance over the maximum permitted weight. (*!!Danger of breaking!!*). Always make sure that there are no living beings or materials below the load that could be injured or damaged!

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